



Reusable Feeding Needle - Barrel Tip - 30 mm/20G/Straight



With luer lock hubs. Ideal for transcardiac perfusions with fixative.

General Info

Product Number	18060-20
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506780



Find 18060-20 on
finescience.com

Technical Info

Shipping Weight	2 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	30 mm		
Gauge	20		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	1.7 mm		
Tipshape	Barrel		
Xsizeweightmc	20-25 g		
Xsizeweightrt	50-75 g		