



Allgemein

Produktnummer	18060-22
Kategorie	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511814



Finde 18060-22 auf
finescience.com

Technische Informationen

Versandgewicht	2 g	Nettogewicht	2 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	25 mm		
Gauge	22		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	1.7 mm		
Tipshape	Barrel		
Xsizeweightmc	15-20 g		