



Reusable Feeding Needle - Barrel Tip - 25 mm/22G/Straight



General Info

Product Number	18060-22
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511814



Find 18060-22 on
finescience.com

Technical Info

Shipping Weight	2 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	25 mm		
Gauge	22		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	1.7 mm		
Tipshape	Barrel		
Xsizeweightmc	15-20 g		