



With luer lock hubs. Ideal for transcardiac perfusions with fixative.



General Info

Product Number	18061-10
Category	Lab Accessories > Feeding Needles
EAN	4251909506797



Find 18061-10 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	5 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	100 mm		
Gauge	16		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	3 mm		
Tipshape	Round		
Xsizeweightrt	150-300 g		
Xsizeweightgp	350-450 g		