



Single Use Feeding Needles - Straight/50mm/18G/2.5mm Tip Diameter



These economical feeding needles feature plastic luer hubs with silicone tips. They are nontoxic and nonpyrogenic. Needles can be easily formed to any desired shape.

General Info

Product Number	18061-18
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506810



Find 18061-18 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Various		
Length	50 mm		
Gauge	18		
Shape	Straight		
Sterile	Yes		
Style	Single use		
Tipdiameter	2.5 mm		
Tipshape	Barrel		