



Reusable Feeding Needle - Ball Tip - 25 mm/22G/Straight/1.25 mm Ø



With luer lock hubs. Ideal for transcardiac perfusions with fixative.

General Info

Product Number	18061-22
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506834



Find 18061-22 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	25 mm		
Gauge	22		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	1.25 mm		
Tipshape	Round		
Xsizeweightmc	15-20 g		