



Reusable Feeding Needle - Ball Tip - 38 mm/18G/Straight/2.25 mm Ø



With luer lock hubs. Ideal for transcardiac perfusions with fixative.

General Info

Product Number	18061-38
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506858



Find 18061-38 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	38 mm		
Gauge	18		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	2.25 mm		
Tipshape	Round		
Xsizeweightmc	20-30 g		
Xsizeweightrt	100-200 g		
Xsizeweightgp	250-300 g		