



Reusable Feeding Needle - Ball Tip - 50 mm/18G/Curved/2.25 mm Ø



With luer lock hubs. Ideal for transcardiac perfusions with fixative.

General Info

Product Number	18061-50
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506865



Find 18061-50 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	50 mm		
Gauge	18		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	2.25 mm		
Tipshape	Round		
Xsizeweightmc	25-30 g		
Xsizeweightrt	100-200 g		
Xsizeweightgp	250-300 g		