



Reusable Feeding Needle - Ball Tip - 75 mm/16G/Curved/3 mm Ø



With luer lock hubs. Ideal for transcardiac perfusions with fixative.

General Info

Product Number	18061-75
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506872



Find 18061-75 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	75 mm		
Gauge	16		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	3 mm		
Tipshape	Round		
Xsizeweightrt	150-300 g		
Xsizeweightgp	350-450 g		