



Kronecker Ligation Aid - 13cm/1mm Hole Diameter



Used for guiding fine suture under and around blood vessels for tying-off purposes, and for ligating peripheral nerves. For item numbers 18062-13 and 18063-13, a groove through the hole helps threading.

General Info

Product Number	18062-13
Category	Needles, Threads & Holders
EAN	4251909506896



Find 18062-13 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	25 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13.5 cm		
Feature	Right handed		
Holediameter	1 mm		
Islefthanded	No		