



Kronecker Ligation Aid - 13cm/1mm Hole Diameter/Left-Handed



Used for guiding fine suture under and around blood vessels for tying-off purposes, and for ligating peripheral nerves. A groove through the hole helps threading.

General Info

Product Number	18063-13
Category	Needles, Threads & Holders
EAN	4251909506902



Find 18063-13 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	27 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13.5 cm		
Holediameter	1 mm		
Islefthanded	Yes		