



Allgemein

Produktnummer	18063-20
Kategorie	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511838



Finde 18063-20 auf
finescience.com

Technische Informationen

Versandgewicht	3 g	Nettogewicht	2 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	38 mm		
Gauge	20		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	1.7 mm		
Tipshape	Barrel		
Xsizeweightmc	20-25 g		
Xsizeweightrt	50-75 g		