



Reusable Feeding Needle - Ball Tip - 38 mm/20G/Curved/2.25 mm Ø



With luer lock hubs. Ideal for transcardiac perfusions with fixative.

General Info

Product Number	18065-20
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909506926



Find 18065-20 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	38 mm		
Gauge	20		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	2.25 mm		
Tipshape	Round		
Xsizeweightmc	20-25 g		
Xsizeweightrt	50-75 g		