



Reusable Feeding Needle - Barrel Tip - 100 mm/16G/Curved



Allgemein

Produktnummer	18067-10
Kategorie	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511845



Finde 18067-10 auf finescience.com

Technische Informationen

Versandgewicht	4 g	Nettogewicht	3 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	100 mm		
Gauge	16		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	2.4 mm		
Tipshape	Barrel		
Xsizeweightrt	150-300 g		
Xsizeweightgp	350-450 g		