



Reusable Feeding Needle - Barrel Tip - 100 mm/16G/Curved



General Info

Product Number	18067-10
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511845



Find 18067-10 on
finescience.com

Technical Info

Shipping Weight	4 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	100 mm		
Gauge	16		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	2.4 mm		
Tipshape	Barrel		
Xsizeweightrt	150-300 g		
Xsizeweightgp	350-450 g		