



Double Pronged Hook - 11cm



The tip consists of two very sharp, nearly horizontally opposed points. By applying gentle downward pressure and twisting slightly in a clockwise fashion, the points will embed themselves horizontally into fine membranes or skin so that it can be lifted upward. This instrument also has been used successfully when separating fine tissue from an underlying area that is not to be cut.

General Info

Product Number	18067-11
Category	Pins, Probes & Hooks > Hooks
EAN	4251909506933



Find 18067-11 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Tipdiameter	1 mm		
Tipshape	Curved		