



Reusable Feeding Needle - Barrel Tip - 38 mm/20G/Curved



General Info

Product Number	18067-20
Category	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511852



Find 18067-20 on
finescience.com

Technical Info

Shipping Weight	2 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	38 mm		
Gauge	20		
Shape	Curved		
Sterile	No		
Style	Reusable		
Tipdiameter	1.7 mm		
Tipshape	Barrel		
Xsizeweightmc	20-25 g		
Xsizeweightrt	50-75 g		