



Horizon™ Ligation Applicator



The Horizon™ System provides a strong, secure clip coupled with an efficient, convenient easy-load cartridge system.

General Info

Product Number	18070-00
Category	Clamps & Clips > Horizon System
EAN	4251909512156



Find 18070-00 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	49 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	15 cm		
Shape	Curved		
Sterile	No		