



Horizon™ Ligation Clips



The preformed chevron shape of the Horizon™ Clip engulfs tissue with precise tip-to-tip closure. Heart-shaped wire is designed to give each clip a firm grip on vessels, while the triangulated cross-section of the clip leg maximizes surface-to-surface contact.

General Info

Product Number	18070-10
Category	Clamps & Clips > Horizon System
EAN	34026704410201



Find 18070-10 on
finescience.com

Technical Info

Shipping Weight	12 g	Net Weight	11 g
Package Size	6 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Autoclavesafe	Yes		
Dimensionstext	1.98 x 2.87 mm		
Sterile	No		