



Vascular Occluder - VO-1.5/3.5mm x 1.5mm



These vascular occluders are ideal for fast and dependable zero-flow baseline determination in blood flow studies. They also provide partial to full vessel occlusion for circulation research studies, and prolonged implantations or constriction of soft organs in acute or chronic applications.

General Info

Product Number	18080-01
Category	Clamps & Clips > Occluders
EAN	4251909506940



Find 18080-01 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Silicone		
Clamppress	Variable		
Jawdims	3.5 x 1.5 mm		
Tipdimensions	A: 1.5 mm - B: < 3.5 mm - C: 1.5 mm		