



S&T Ball Chain Retraction System - Ball Chain Elevators



The Base Plate is made of matte stainless steel to reduce glare. The perimeter slots will secure the Ball Chain Retractors, offering retraction in every direction (lateral and upwards) with the help of the Ball Chain Elevators. The Ball Chain Retractors can also be used to hold the animal firmly in position. There are four feet which lift the Base Plate by 1.5cm to facilitate the locking of the Ball Chain Retractors. The Wide Hooks can be attached to the Ball Chain Retractors if desired.

This product is discontinued and has limited availability.

General Info

Product Number	18100-02
Category	Retractors > Retraction Systems
EAN	4251909509552



Find 18100-02 on finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	63 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Width	3 cm		
Height	8.5 cm		
Dimensionstext	8.5cm High x 3cm Wide		