



Absorbent Swabs - Mounted/Sterile



These cellulose surgical sponges absorb fluids twice as fast as the best cotton sponges and have three times their fluid-holding capacity. They are made from strong, soft, lint-free, biocompatible-proven material which is non-irritating to tissue and leaves no fibrous residue which might cause wound contamination. They are sterilized by radiation and come in clear, sealed packages.

General Info

Product Number	18105-01
Category	Lab Accessories Most Popular > Lab Accessories
EAN	4251909506988



Find 18105-01 on
finescience.com

Technical Info

Shipping Weight	225 g	Net Weight	140 g
Package Size	200 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Cellulose with plastic handle		
Length	7 cm		
Mounted	Yes		
Sterile	Yes		