



Absorbent Swabs - Unmounted / Sterile



These cellulose surgical sponges absorb fluids twice as fast as the best cotton sponges and have three times their fluid-holding capacity. They are made from strong, soft, lint-free, biocompatible-proven material which is non-irritating to tissue and leaves no fibrous residue which might cause wound contamination.

General Info

Product Number	18105-04
Category	Lab Accessories
EAN	4251909511760



Find 18105-04 on
finescience.com

Technical Info

Shipping Weight	153 g	Net Weight	10 g
Package Size	200 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Cellulose		
Length	20 mm		
Mounted	No		
Sterile	Yes		