



Vessel Dilating Forceps - Angled 10°/12cm/0.3mm Tip Diameter



Highly modified forceps with elongated, parallel tips which are highly polished. Designed for controlled intraluminal vessel dilation.

This product is discontinued and has limited availability.

General Info

Product Number	18133-12
Category	Forceps & Tweezers > Smooth
EAN	4251909509590



Find 18133-12 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	12.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdiameter	0.3 mm		
Tips	For Dilating		
Tipshape	Angled 10°		