



Mini Bulldog Artery Clip - Straight/30mm



This bulldog artery spring clamp provides minimum distortion during MRI procedures.

General Info

Product Number	18144-30
Category	Clamps & Clips > Clamps
EAN	4251909507039



Find 18144-30 on finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Clampinglength	9 mm		
Clamppress	Strong		
Length	30 mm		
Grampress	230 g		
Jawdims	9 x 1 mm		
Serrated	Yes		
Shape	Straight		