



Pierse Fixation Forceps - Straight/13cm



A 0.3mm tissue retention hole provides an atraumatic hold for sclera, vascular sheath and epineurium.

General Info

Product Number	18155-13
Category	Forceps & Tweezers > Smooth
EAN	4251909507046



Find 18155-13 on
finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	1 x 0.5 mm		
Tipwidth	1.5 mm		
Tips	0.3 mm Hole		
Tipshape	Straight		