



This Small Animal Retraction System for animals from 200 g to 5 kg simplifies set-up procedure and provides the examiner with an excellent overview of the operative site. It consists of three components:

1. Base Plate (Nos. 18200-03 & 18200-04)
2. Fixators (Nos. 18200-01 & 18200-02)
3. Retractors (Nos. 18200-05 through 18200-12)



The Base Plates are made from magnetic stainless steel and are available in two sizes. The surface is nonreflective and both plates have a window in the center which allows the animal contact with body temperature control systems.

The window sizes are 10 x 13cm for the small base plate and 15 x 25cm for the large.

The Fixators are movable anchors made from anodized aluminum with a surgical grade rare earth magnet in the base for attachment anywhere on the Base Plate.

The magnets do not interfere with lab electronics and can be used in close proximity (3cm) to ferro-magnetic items.

They incorporate pushbutton-activated spring locks that grip either the Wire or Elastomer Retractors.

The Retractors are interchangeable Retractor Tips (Nos. 18200-08 through 18200-12) which attach to either Retractor Wire (Nos. 18200-05 & 18200-06) for static retraction or to Elastomer (No. 18200-07) for dynamic retraction.

The Retractor Tips are light and atraumatic to delicate tissue or membrane.

All components of the system are autoclavable, with the exception of the elastomers (No. 18200-07)

Components can be purchased individually or in kits.

Kit for animals from 200g up to 5kg (No. 18200-50) includes:

- 1 Large Base Plate (No. 18200-04) 2 Short Fixators (No. 18200-01)
- 2 Tall Fixators (No. 18200-02)
- 1 Elastomer, 2m roll (No. 18200-07)
- 1 package of Blunt Retractors, 2.5mm wide (No. 18200-10)
- 1 package of Blunt Retractors, 5mm wide (No. 18200-11)

General Info

Product Number	18200-50
Category	Retractors > Retraction Systems Most Popular > Retractors
EAN	4251909507183



Find 18200-50 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	775 g	Net Weight	770 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Feature	Kit for animals from 200 g to 5 kg		