



### Allgemein

Produktnummer	18320-11
Kategorie	Clamps & Clips > Clamps
EAN	4251909507190



Finde 18320-11 auf  
[finescience.com](https://www.finescience.com)

### Technische Informationen

Versandgewicht	5 g	Nettogewicht	1 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Clampinglength	13 mm		
Clamppress	Light		
Length	33 mm		
Grampress	15 g		
Jawdims	13 x 1.5 mm		
Serrated	Yes		
Shape	Straight		