

F · S · T

Moria Vessel Clamp - Straight/33mm



General Info

Product Number	18320-11
Category	Clamps & Clips > Clamps
EAN	4251909507190



Find 18320-11 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	13 mm		
Clamppress	Light		
Length	33 mm		
Grampress	15 g		
Jawdims	13 x 1.5 mm		
Serrated	Yes		
Shape	Straight		