



Moria De Wecker Clamp - Curved/50mm



Allgemein

Produktnummer	18374-44
Kategorie	Clamps & Clips > Clamps
EAN	4251909507213



Finde 18374-44 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	3 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	15 mm		
Clamppress	Strong		
Length	50 mm		
Jawdims	16 x 2.5 mm		
Serrated	No		
Shape	Curved		