



Vessel Cannulation Forceps - 0.35mm OD/11.5cm



These vessel cannulation forceps are useful whenever a fine plastic tube has to be introduced into a small blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is directed exactly into the vessel opening without damage to the vessel from unwanted movement.

General Info

Product Number	18403-11
Category	Forceps & Tweezers > Cannulation
EAN	4251909507220



Find 18403-11 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Tipshape	Angled		
Tubingsize	0.35 mm OD		