



Vessel Cannulation Forceps - 0.6mm OD/13.5cm



These vessel cannulation forceps are useful whenever a fine plastic tube has to be introduced into a small blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is directed exactly into the vessel opening without damage to the vessel from unwanted movement.

General Info

Product Number	18406-13
Category	Forceps & Tweezers > Cannulation
EAN	4251909507244



Find 18406-13 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	55 g	Net Weight	32 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13.5 cm		
Tipshape	Angled		
Tubingsize	0.6 mm OD		