



Micro Clips - Titanium/Straight/19mm



Lightweight with strong clamping pressure. Delicate atraumatic serrations.

General Info

Product Number	18555-01
Category	Clamps & Clips > Serrefines
EAN	4251909507251



Find 18555-01 on finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Clamppress	Light		
Length	19 mm		
Grampress	60 g		
Jawdims	10 x 2 mm		
Serrated	Yes		
Shape	Straight		