



Micro Clips - Titanium/Straight/15mm



Lightweight with strong clamping pressure. Delicate atraumatic serrations.

General Info

Product Number	18555-03
Category	Clamps & Clips > Serrefines
EAN	4251909507275



Find 18555-03 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Clamppress	Medium		
Length	15 mm		
Grampress	100 g		
Jawdims	6 x 1 mm		
Serrated	Yes		
Shape	Straight		