



Micro Clips - Tanium/Straight/13mm



Lightweight with strong clamping pressure. Delicate atraumatic serrations.

General Info

Product Number	18555-04
Category	Clamps & Clips > Serrefines
EAN	4251909507282



Find 18555-04 on finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Clamppress	Strong		
Length	13 mm		
Grampress	125 g		
Jawdims	4 x 0.75 mm		
Serrated	Yes		
Shape	Straight		