



Multipurpose Forceps - Straight/Suture Tying Platform/14.5cm/0.1mm
Tip Diameter



Balanced instruments are designed to put the center of gravity between the thumb and index finger. This is accomplished by adding a slight counter weight to the proximal end of the instruments. The design enhances tactile feel, aiding in better control and precision. These multipurpose forceps can be used in different surgical techniques. The tips are configured to enable their use as dilators at the front but also needle holders and tying forceps further in the back. The platform at the tips close parallel ensuring a reliable grip.

General Info

Product Number	18601-15
Category	Forceps & Tweezers > Smooth
EAN	4251909507299



Find 18601-15 on
finescience.com

Technical Info

Shipping Weight	65 g	Net Weight	30 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdiameter	0.1 mm		
Tips	For Dilating and Suture Tying		
Tipshape	Straight		