



Multipurpose Forceps - Straight/Suture Tying Platform/15cm/0.3mm Tip Diameter



Balanced instruments are designed to put the center of gravity between the thumb and index finger. This is accomplished by adding a slight counter weight to the proximal end of the instruments. The design enhances tactile feel, aiding in better control and precision. These multipurpose forceps can be used in different surgical techniques. The tips are configured to enable their use as dilators at the front but also needle holders and tying forceps further in the back. The platform at the tips close parallel ensuring a reliable grip.

General Info

|                |                             |
|----------------|-----------------------------|
| Product Number | 18602-15                    |
| Category       | Forceps & Tweezers > Smooth |
| EAN            | 4251909507305               |



Find 18602-15 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 65 g                            | Net Weight      | 32 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 15 cm                           |                 |      |
| Tipdiameter            | 0.3 mm                          |                 |      |
| Tips                   | For Dilating and Suture Tying   |                 |      |
| Tipshape               | Straight                        |                 |      |