



Single Use Feeding Needles (PTFE) - 85 mm/17/2.4 mm



General Info

Product Number	18961-85
Category	Lab Accessories > Feeding Needles
EAN	4251909510329



Find 18961-85 on
finescience.com

Technical Info

Shipping Weight	11 g	Net Weight	10 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Polytetrafluoroethylene (PTFE)		
Length	85 mm		
Gauge	17		
Shape	Flexible		
Sterile	Yes		
Style	Single use		
Tipdiameter	2.4 mm		
Tipshape	Barrel		