



Gold Wire Contacts - Female

Matching pin and socket. Contacts are crimp terminated, but also can be used for solder connections.



General Info

Product Number	19003-01
Category	Lab Accessories
EAN	4251909507336



Find 19003-01 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	50 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Diameter	1.2 mm		
Length	1 cm		
Lengthwhenconnected	16 mm		
Style	Female		
Wireopening	0.63 mm		