



Burrs for Micro Drill - 0.9mm Tip Diameter/Carbon Steel



MOST
POPULAR



NOT
AUTOCLAVABLE



General Info

Product Number	19007-09
Category	Lab Accessories > Micro Drill Accessories Most Popular > Lab Accessories
EAN	4251909507367



Find 19007-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	8 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Length	44 mm		
Shaftdiameter	2.3 mm		
Tipdiameter	0.9 mm		