



Burrs for Micro Drill - 1.4mm Tip Diameter/Carbon Steel



MOST
POPULAR



NOT
AUTOCLAVABLE



Allgemein

Produktnummer	19007-14
Kategorie	Lab Accessories > Micro Drill Accessories Most Popular > Lab Accessories
EAN	4251909507374



Finde 19007-14 auf
finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	9 g
Paketgröße	10 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Length	44 mm		
Shaftdiameter	2.3 mm		
Tipdiameter	1.4 mm		