



Burrs for Micro Drill - 0.9mm Tip Diameter/Stainless Steel



Allgemein

Produktnummer	19008-09
Kategorie	Lab Accessories > Micro Drill Accessories
EAN	4251909507411



Finde 19008-09 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	8 g
Paketgröße	10 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	44 mm		
Shaftdiameter	2.3 mm		
Tipdiameter	0.9 mm		