



Burrs for Micro Drill - Tungsten Carbide/Yes/0.5 mm



General Info

Product Number	19009-05
Category	Lab Accessories > Micro Drill Accessories
EAN	4251909507442



Find 19009-05 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Tungsten Carbide		
Autoclavesafe	Yes		
Length	44 mm		
Shaftdiameter	2.3 mm		
Tipdiameter	0.5 mm		