



Burrs for Micro Drill - Tungsten Carbide/Yes/0.9 mm



General Info

Product Number	19009-09
Category	Lab Accessories > Micro Drill Accessories
EAN	4251909507466



Find 19009-09 on [finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Tungsten Carbide		
Autoclavesafe	Yes		
Length	44 mm		
Shaftdiameter	2.3 mm		
Tipdiameter	0.9 mm		