



Self-Tapping Bone Screws - 4.8mm/100 screws/1.19mm Shaft Diameter

These screws have been found to hold better than ordinary machine screws. They have a fillister flat head and are ideal for anchoring dental cement to bone. Use our 0.9mm burs No. 19007-09 for pre-drilling the 1.19mm No. 19010-00 screws.



General Info

Product Number	19010-00
Category	Lab Accessories > Screws
EAN	4251909507473



Find 19010-00 on finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	4.8 mm		
Feature	100 screws		
Shaftdiameter	1.19 mm		