



Self-Tapping Bone Screws - 4mm/100 screws/0.85mm Shaft Diameter



These screws have been found to hold better than ordinary machine screws. They have a fillister flat head and are ideal for anchoring dental cement to bone. Use our 0.7mm burs No. 19007-07 for pre-drilling the 0.86mm (No. 19010-10) screws.

General Info

Product Number	19010-10
Category	Lab Accessories > Screws
EAN	4251909507480



Find 19010-10 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	4 mm		
Feature	100 screws		
Shaftdiameter	0.86 mm		