



Self-Tapping Bone Screws - 4 mm/10 screws/0.86 mm Shaft Diameter



Allgemein

Produktnummer	19010-11
Kategorie	Lab Accessories > Screws
EAN	4251909510343



Finde 19010-11 auf finescience.com

Technische Informationen

Versandgewicht	6 g	Nettogewicht	5 g
Paketgröße	10 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	4 mm		
Feature	10 screws		
Shaftdiameter	0.86 mm		