



Moria Instrument Case with Silicone Mat - 18 x 9 x 3cm

This instrument case comes with a spiked silicone mat for secure storage of surgical instruments. The case and silicone insert are autoclavable but should not be used for sterilization of surgical instruments.



General Info

Product Number	20311-21
Category	Instrument Care > Storage
EAN	4251909507497



Find 20311-21 on
finescience.com

Technical Info

Shipping Weight	220 g	Net Weight	210 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Aluminum		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Color	Silver		
Dimensionstext	18 x 9 x 3 cm		