



Colored Aluminum Instrument Case - 18 x 9 x 3 cm/Yellow



This instrument case includes a spiked silicone mat for secure storage of surgical instruments. It is available in five colors, though slight variations may occur between production batches. Both the case and the silicone insert are autoclavable, but they are not intended for use in the sterilization of surgical instruments.

General Info

Product Number	20330-03
Category	Instrument Care > Storage Most Popular > Instrument Care
EAN	4251909511517



Find 20330-03 on
finescience.com

Technical Info

Shipping Weight	207 g	Net Weight	207 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Aluminum		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Color	Yellow		
Length	18 cm		
Width	9 cm		
Height	3 cm		