




Plastic Sterilization Case - 27 x 16 x 4cm



This instrument case comes with a spiked silicone mat for secure storage of surgical instruments. Perforated top is clear for visual identification of instruments. Suitable for autoclave or ethylene oxide sterilization.

General Info

Product Number	20810-02
Category	Instrument Care > Sterilization
EAN	4251909507541



Find 20810-02 on
finescience.com

Technical Info

Shipping Weight	655 g	Net Weight	644 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Dimensionstext	27 x 16 x 4 cm		